

Online Graph: 219

FCU 2210-4 U3.29

HYDAC FCU 2210/4

Start: 16.04.2010 10:10

Averaging interval: 0 min


Test volume : 100 ml

ISO > 4µm : _____

ISO > 6µm : _____

ISO > 14µm :

0 10 12 14 16 18 20 22



Last Measurement:

size	counts per 100ml			
µm	diff	cum	SAE	ISO
4	2831	3887	03	12
6	723	656	02	10
14	61	133	02	08
21	72	72	03	